

## HT-3F(5)

Electric Operating Table  
(Five Functions)



- This table is suitable for general, gynecology, obstetric, proctology, urology and orthopaedics operation.
- Radiolucent table top and longitudinal moving 300mm function are very convenient for X-ray examination and C-arm photograph on different parts of body.
- It adopts Jiechang brand motors, smooth and silent operation. The height, back section, trendelenburg, lateral and moving 300mm adjustments are all controlled by motors.
- The head and leg section are manual adjustment.
- Made of full premium 304 stainless steel. Mobile stainless steel base on castors with central brake.
- All sections fit with mattress, detachable for easy cleaning.

### Accessories:

Mattress : 1pc  
Power line: 1pc  
Arm supports : 2pcs  
Leg supports : 2pcs  
Hand controller: 1pc  
Waist supports : 2pcs  
Shoulder supports : 2pcs  
Controller on the column: 1pc  
Anesthetic screen frame : 1pc

**Optional :** Two types of base and two types of mattress

Size	2050*500mm	Back Section (electric)	70°/30°
Height Adjustment (electric )	700-1000mm	Head Board (manual)	45°/90°
Trendelenburg/Reverse Trendelenburg (electric )	25°	Kidney Bridge	0-120mm
Left/Right tilt (electric )	20°	Leg Board (manual )	15°/90°
		Table Top Longi tudinal Movi ng(electri c)	300mm